



Sodium Hyaluronate 1.8% w/v (Can also be supplied in 1.0%)

Salient Features :

- Supplied Sterile in a 1ml Pre Filled Syringe with a 27G cannula
- This is a sterile, non-pyrogenic, viscoelastic, preparation, of sodium hyaluronate. It contains 18mg/ml of sodium hyaluronate dissolved in physiological sodium chloride phosphate buffer. (Other concentrations like 1.0% also available)
- Sodium Hyaluronate is derived from bacterial fermentation. Its viscosity is responsible for maintenance of space in the anterior chamber and allows the manipulation of tissue during the surgery without escape through the incision.
- The 1.8% solution is transparent . It does not interfere with epithelization and normal wound healing.
- Uses : It is indicated for use as a surgical aid in cataract extraction (intra and extra capsular extraction), IOL implant, corneal transplant, glaucoma filtration and retinal detachment surgery.
- Contraindication : At present there are no known contraindications of the use of this sodium hyaluronate.
- Cost economical.

Another Quality Product from :



Action Medical Mktg Pvt Ltd.

(An ISO 9001:2008 & ISO 13485:2003 Certified Company)

Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road

Taluka Shirur, Dist Pune Pin 412216, INDIA

Tel : 0091-2137-253617 / 18. Fax : 0091-2137-253620.

E-mail : action@actionmedihealth.com

Website : www.actionmedihealth.com