

## EYE DRAPE - OPHTHALMIC

Model : **DISP-No 01**



**Eye Drape – OPHTHALMIC**  
**Brand Action / PrecisionCut**

Convenient to use, hygienic & sterile, this range of disposable ophthalmic supplies that we offer. These ACTION ophthalmic supplies are meticulously designed with attention to doctor's & patient's need for convenient use. The usage of qualitative raw materials used assures high performance, safety and affordability.

These are plastic drapes, folded and supplied sterile. Meant for Ophthalmic surgeries.

Applications: Used as a drape for providing a precise working area for ophthalmic surgeries.

Best Suitable: Ophthalmic procedures.

Product's USP:

- Fits both Left or Right eye.
- Soft Plastic drape-able material.
- **Specs :**
  - Towel Area 60x60cms or 46x60cms or 40x40cms.
  - Adhesive Area : 6x6cms
  - With Drain pouch.
  - Supplied Sterile in individual packs.

Another Quality Product from :



**Action Medical Mktg Pvt Ltd.**

(An ISO 9001:2008 & ISO 13485:2003 Certified Company )

Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road

Taluka Shirur, Dist Pune Pin 412216, INDIA

Tel : 0091-2137-253617 / 18. Fax : 0091-2137-253620.

E-mail : [action@actionmedihealth.com](mailto:action@actionmedihealth.com)

Website : [www.actionmedihealth.com](http://www.actionmedihealth.com)