


ACTION STERILIZATION TRAYS

For the instruments you value most



Available in Acrylic & Stainless Steel in various sizes

Acrylic Small Size (IM31103)	: (15 x 7 x 2)cm	US\$16
Acrylic Medium Size (IM31105)	: (20 x 10 x 3.5)cm	US\$40
Acrylic Large Size (IM31109)	: (28 x 15 x 5) cm	US\$140

Stainless Steel Trays

Small (IM3101)	: (21x 12x2) cm	US\$24
Medium (IM3102)	: (25x21x 2) cm	US\$56
Large (IM3104)	: (42x25x 2) cm	US\$180

ACTION Sterilization Trays are designed to address specific requirements of ophthalmologists and other surgeons who are using delicate and precious instruments for microsurgeries.

All instruments are not equal. Some are special. **ACTION Sterilization Tray** is meant for those instruments that you value the most, your diamond knives, Mcpherson forceps or corneal scissors to name just a few. Compact, Light and easy to carry, **ACTION Sterilization Tray** enables you to easily monitor your precious possessions before, during and after sterilization.

Remember, your surgery depends on them.

SALIENT FEATURES (Acrylic) :

- Made of superior quality plastic thus eliminates the possibility of wearing off the sharpness of delicate instruments.
- Can withstand typical repetitive wet autoclaving process cycle with maximum temperature of 121 °C and a pressure of 15 psi.
- Transparent lid enables you to identify the instruments without opening the box and thus eliminating possibility of compromising sterility anytime prior to surgery.
- Snap-in lid allows trays to be accommodated in the available space within the sterilizer without a possibility of inadvertent opening.
- Silicon Mat holds instruments in their place without rubbing against each other or to the container walls.



Action Medical Pvt Ltd.

(An ISO 9001:2008 & ISO 13485:2003 Certified Company)

Gat No 528, Koregaon-Bhima, Near SBI, Pune-Nagar Road

Taluka Shirur, Dist Pune Pin 412216, INDIA

Tel : 0091-2137-253617 / 18. Fax : 0091-22137-253620.

E-mail : action@actionmedihealth.com

Website : www.actionmedihealth.com